

AbsoluteBMC™
Bone Marrow Concentrating System
GenesisCS Component Concentrating System
Date: May 2024
Instruction for Use

ATTENTION OPERATING SURGEON

NOTE: DEVICE IS FOR SINGLE USE ONLY. Discard the entire disposable system after one use, using an acceptable disposal method for products potentially contaminated with blood.

DESCRIPTION

1. The AbsoluteBMC™ GenesisCS Component Concentrating System is manufactured by EmCyte Corporation. The kit prepares platelet poor plasma and platelet concentrate from a small sample of blood and a cell concentrate from bone marrow at the point of care. The system contains syringes, a bone marrow needle accessory, and the concentrating device accessories.

MATERIALS

2. The materials used are syringes, needles, tubing, connectors, and concentrating devices. The materials consist of medical grade polymers, elastomers and stainless steel that are suitable for use in medical devices. All components in this system are packaged, labeled and sterilized as indicated by the manufacturer's labeling. All components in this system are latex-free.

INDICATIONS FOR USE STATEMENTS

3. The GenesisCS Component Concentrating System is intended to be used in a clinical laboratory or intraoperatively at the point of care for the safe and rapid preparation of platelet poor plasma and platelet concentrate from a small sample of blood and for a preparation of a cell concentrate from bone marrow. The safety and effectiveness of this device for in vivo indications for use has not been established.
4. The safety and effectiveness of this device for in vivo indications for use, such as bone healing and hemostasis, have not been established.
5. The PRP and BMC prepared by this device has not been evaluated for any clinical indications.
6. The PRP and BMC prepared by this device is NOT indicated for delivery to the patient's circulatory system.

USER POPULATION

7. The intended user population is medical professionals who are licensed
8. or certified in clinical practice. The operational context of the device
9. requires users to be trained on aseptic technique and understand blood
10. components. The surgeon is to be thoroughly familiar with the
11. equipment and the surgical procedure prior to using this device.

DEVICE USE ENVIRONMENT

12. The device is intended to be used in a health care setting such as a surgery room, clinic or outpatient care center.

WARNING AND PRECAUTIONS

13. Use proper safety precautions to guard against needle sticks.
14. Follow manufacturer instructions when using centrifuge. Use only EmCyte provided general purpose centrifuge. Outcomes using centrifuges from other manufacturers are unknown.
15. When using the bone marrow aspiration needle, follow manufacturer's instruction for use.
16. Do not use sterile components of this system if package is opened or damaged.
17. Single use device. Do not reuse. Do not attempt to clean or re-sterilize this product.
18. Do not use after expiration date.
19. Use prepared BMC, PPP or PRP within 4 hours after drawing blood or bone marrow aspirate according to current AABB guidelines.
20. BMC prepared from bone marrow may contain higher levels of plasma free hemoglobin than PRP prepared from whole blood.

POSSIBLE RISKS

21. The patient is to be made aware of the general risks associated with bone marrow aspiration. These risks include, but are not limited to: hemorrhage, seroma formation, infection, and/or persistent pain at the site of aspiration.
22. Reuse may be a potential biohazard

POSSIBLE ADVERSE EFFECTS

23. Damage to blood vessels, hematoma, delayed wound healing and/or infection is associated with blood draw, bone marrow harvest and/or surgical procedure.
24. Temporary or permanent nerve damage that may result in pain or numbness is associated with blood draw, and/or surgical procedure.
25. Early or late postoperative infection is associated with surgical procedure.
26. Pain associated with site of bone marrow harvest.

STERILITY

27. The AbsoluteBMC™ GenesisCS Component Concentrating System kits are sterilized by ETO exposure. Do not use any component from an opened or damaged package. Do not resterilize. Discard if kit packaging is damaged or open.

INSTRUCTIONS FOR USE FOR 60mL SYSTEM

PREPARATION PROTOCOL

NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.

28. **BONE MARROW ASPIRATE DRAW:**
Attach the sterile filter needle onto the 60mL VACLOK syringe. Draw 15mL of Heparin Anticoagulant (1000 units/mL) into the 60mL syringe. Remove the filter needle from the syringe. Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it. Attach the heparin syringe to the OUT port of the bone marrow filter. Inject Heparin filling the filter to prime and then aspirate the heparin back into the syringe. Then discard heparin leaving 5mL in the VACLOK syringe. Slowly draw 55mL of bone marrow aspirate (BMA) from the patient, filling the syringe to 60mL. Follow the bone marrow needle manufacturer's package insert to obtain bone marrow aspirate. Gently, but thoroughly mix the BMA and heparin upon collection to prevent coagulation.
29. **FILTER:** Connect a 60mL syringe to the OUT port of the bone marrow filter. Connect the VACLOK bone marrow syringe to the IN port of the filter. Inject the BMA through the filter into the 60mL syringe. Once completed, clear the remaining bone marrow in the filter by aspirating it into the 60mL syringe.

CONCENTRATING PROTOCOL

30. **LOAD:** Remove and discard the red vented cap from the needle-less port of the Concentrating Device. Slowly add the filtered and anticoagulated BMA through the port of the Concentrating Device.
31. **BALANCE:** Make sure the counterbalance device contains the same amount of volume as the Concentrating Device. Then place them directly opposite to each other in the centrifuge rotor buckets.
32. **CENTRIFUGATION:**
 - a. Sapphire Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - b. Platinum Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - c. Executive Series Centrifuge: Close the lid and set to 10 minutes and 4.4 x 1000 RPM (4400 RPM).
 - d. Press the start button. Once the centrifuge stops, remove the Concentrating Device.
33. **BMC EXTRACTION:** Attach the sterile 60mL syringe to the needle-less port and aspirate the plasma until the bone marrow RBC interface reaches the 7mL volume marker, then stop aspirating. Then attach the 12mL syringe and aspirate 7mL of BMC. Remove sterile syringe and apply a sterile cap.

INSTRUCTIONS FOR USE FOR 120mL SYSTEM

INSTRUCTIONS FOR USE FOR 30mL SYSTEM

PREPARATION PROTOCOL

NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.

34. **BONE MARROW ASPIRATE DRAW:**
Attach the sterile filter needle onto the 30mL VACLOK syringe. Draw 15mL of Heparin Anticoagulant (1000 units/mL) into the 30mL syringe. Remove the filter needle from the syringe. Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it. Attach the heparin syringe to the OUT port of the bone marrow filter. Inject Heparin filling the filter to prime and then aspirate the heparin back into the syringe. Then discard heparin leaving 3mL in the VACLOK syringe. Slowly draw 27mL of bone marrow aspirate (BMA) from the patient, filling the syringe to 30mL. Follow the bone marrow needle manufacturer's package insert to obtain bone marrow aspirate. Gently, but thoroughly mix the BMA and heparin upon collection to prevent coagulation.
35. **FILTER:** Connect a 30mL syringe to the OUT port of the bone marrow filter. Connect the VACLOK bone marrow syringe to the IN port of the filter. Inject the BMA through the filter into the 30mL syringe. Once completed, clear the remaining bone marrow in the filter by aspirating it into the 30mL syringe.

CONCENTRATING PROTOCOL

36. **LOAD:** Remove and discard the red vented cap from the needle-less port of the Concentrating Device. Slowly add the filtered and anticoagulated BMA through the port of the Concentrating Device.
37. **BALANCE:** Make sure the counterbalance device contains the same amount of volume as the Concentrating Device. Then place them directly opposite to each other in the centrifuge rotor buckets.
38. **CENTRIFUGATION:**
 - a. Sapphire Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - b. Platinum Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - c. Executive Series Centrifuge: Close the lid and set to 10 minutes and 4.4 x 1000 RPM (4400 RPM).
 - d. Press the start button. Once the centrifuge stops, remove the Concentrating Device.
39. **BMC EXTRACTION:** Attach the sterile 30mL syringe to the needle-less port and aspirate the plasma until the bone marrow RBC interface reaches the 3mL volume marker, then stop aspirating. Then attach the 12mL syringe and aspirate 3-4mL of BMC. Remove sterile syringe and apply a sterile cap.

PREPARATION PROTOCOL













NOTE: Use standard sterile aseptic technique throughout the following procedure. Always swab needle-less ports with alcohol before and after accessing.

40. **BONE MARROW ASPIRATE DRAW:**
Attach the sterile filter needles onto the 60mL VACLOK syringes. Draw 15mL of Heparin Anticoagulant (1000 units/mL) into each of the 60mL syringes. Remove the filter needle from the syringe. Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it. Attach the heparin syringe to the OUT port of the bone marrow filter. Inject 5mL into the filter to prime and then aspirate the heparin back into the syringe. Then discard heparin leaving 5mL in the VACLOK syringe. Then attach the sterile filter needle onto the second VACLOK syringe and draw 5mL of Heparin Anticoagulant (1000 units/mL) into it. For each VACLOK syringe, slowly draw 55mL of bone marrow aspirate (BMA) from the patient, filling each syringe to 60mL. Follow the bone marrow needle manufacturer's package insert to obtain bone marrow aspirate. Gently, but thoroughly mix the BMA and heparin upon collection to prevent coagulation.
41. **FILTER:** For each VACLOK syringe do the following steps. Connect a 60mL syringe to the OUT port of the bone marrow filter. Connect the VACLOK bone marrow syringe to the IN port of the filter. Inject the BMA through the filter into the 60mL syringe. Once completed, clear the remaining bone marrow in the filter by aspirating it into the 60mL syringe.

CONCENTRATING PROTOCOL

42. **LOAD:** For each Concentrating Device do the following steps. Remove and discard the red vented cap from the needle-less port of the Concentrating Device. Slowly add the filtered and anticoagulated BMA through the port of the Concentrating Device.
43. **BALANCE:** Make sure each Concentrating device contains the same amount of volume. Then place them directly opposite to each other in the centrifuge rotor buckets.
44. **CENTRIFUGATION:**
 - a. Sapphire Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - b. Platinum Series Centrifuge: Close the lid and set to ABSOLUTEPRP/BMC.
 - c. Executive Series Centrifuge: Close the lid and set to 10 minutes and 4.4 x 1000 RPM (4400 RPM).
 - d. Press the start button. Once the centrifuge stops, remove the Concentrating Device.
45. **BMC EXTRACTION:** For each device, attach the sterile 60mL syringe to the needle-less port and aspirate the plasma until the bone marrow RBC interface reaches the 7mL volume marker, then stop aspirating. Then attach the 12mL syringe and aspirate 7mL of BMC. Collect a total of 14mL of BMC. Remove sterile syringe and apply a sterile cap.





Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

 Do not use if package is damaged	 Single use only	 EmCyte Corporation 4331 Veronica S. Shoemaker Blvd. Fort Myers, FL 33916 Phone: 239-481-7725	 Store in a cool place
 STERILE EO	 Do not re-sterilize	 Attention, read instruction for use	 Store in a dry place
 Consult instruction for use	 Rx Only Prescription Use	 Non-pyrogenic	 MD Medical Device






GSBMA-30: IFU ILLUSTRATION

**NOTICES: PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE.
ALWAYS SWAB SELF-SEALING PORT WITH STERILE ALCOHOL PRIOR TO ACCESSING WITH A STERILE SYRINGE.**

PREPARATION PROTOCOL

<p>STEP 1:</p>  <p>Attach the sterile filter needle onto the VACLOK 30mL syringe. Then draw 15mL of Heparin Anticoagulant</p> <p>(1000 units/mL)</p>	<p>STEP 2:</p>  <p>Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it</p>	<p>STEP 3:</p>  <p>Attach to the OUT port of the bone marrow filter. Fill to prime and then aspirate back into the syringe</p>	<p>STEP 4:</p>  <p>Then discard the residual heparin leaving 3mL in the VACLOK syringe</p>
---	---	---	--

CONCENTRATING PROTOCOL

<p>STEP 1:</p>  <p>Slowly draw 27mL of bone marrow aspirate, filling the syringe to 30mL. Mix the BMA and heparin upon collection to prevent coagulation.</p>	<p>STEP 2:</p>  <p>Start End</p> <p>Connect a 30mL syringe to the OUT port of the bone marrow filter. Connect the VACLOK bone marrow syringe to the IN port of the filter. Inject the BMA through the filter into the 30mL syringe.</p>	<p>STEP 3:</p>  <p>Inject anticoagulated & filtered BMA into the Concentrating Device</p>	<p>STEP 4:</p>  <p>Counterbalance with the same volume as the Concentrating Device.</p>  <p>Then place them directly opposite to each other in the centrifuge rotor buckets.</p>
---	---	---	---


STEP 5:

Sapphire Series
Centrifuge:
AbsolutePRP/BMC

Platinum Series
Centrifuge:
AbsolutePRP/BMC


Executive Series
Centrifuge:
10 minutes
4400 RPM

STEP 6:



Attach the 30mL syringe and aspirate the plasma until the bone marrow RBC interface reaches the 3mL volume marker, then stop aspirating

STEP 7:







Attach the 12mL syringe and aspirate 3-4mL BMC
(Rotate the syringe to re-suspend the BMC buffycoat into the plasma)






GSBMA-60: IFU ILLUSTRATION

**NOTICES: PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE.
ALWAYS SWAB SELF-SEALING PORT WITH STERILE ALCOHOL PRIOR TO ACCESSING WITH A STERILE SYRINGE.**

PREPARATION PROTOCOL:

<p>STEP 1:</p>  <p>Attach the sterile filter needle onto the VACLOK 60mL syringe. Then draw 15mL of Heparin Anticoagulant (1000 units/mL)</p>	<p>STEP 2:</p>  <p>Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it</p>	<p>STEP 3:</p>  <p>Attach to the OUT port of the bone marrow filter. Fill to prime and then aspirate back into the syringe</p>	<p>STEP 4:</p>  <p>Then discard the residual heparin leaving 5mL in the VACLOK syringe</p>
--	---	---	--

CONCENTRATING PROTOCOL

<p>STEP 1:</p>  <p>Slowly draw 27mL of bone marrow aspirate, filling the syringe to 30mL. Mix the BMA and heparin upon collection to prevent coagulation.</p>	<p>STEP 2:</p>  <p>Start End</p> <p>Connect a 30mL syringe to the OUT port of the bone marrow filter. Connect the VACLOK bone marrow syringe to the IN port of the filter. Inject the BMA through the filter into the 30mL syringe.</p>	<p>STEP 3:</p>  <p>Inject anticoagulated & filtered BMA into the Concentrating Device</p>	<p>STEP 4:</p>  <p>Counterbalance with the same volume as the Concentrating Device.</p>  <p>Then place them directly opposite to each other in the centrifuge rotor buckets.</p>
---	---	---	---


STEP 5:

Sapphire Series
Centrifuge:
AbsolutePRP/BMC

Platinum Series
Centrifuge:
AbsolutePRP/BMC


Executive Series
Centrifuge:
10 minutes
4400 RPM

STEP 6:



Attach the 60mL syringe and aspirate the plasma until the bone marrow RBC interface reaches the 7mL volume marker, then stop aspirating

STEP 7:







Attach the 12mL syringe and aspirate 7mL BMC
(Rotate the syringe to re-suspend the BMC buffycoat into the plasma)






GSBMA-120: IFU ILLUSTRATION

**NOTICES: PLEASE DISCARD RED VENTED CAP FROM CONCENTRATING DEVICE BEFORE USE.
ALWAYS SWAB SELF-SEALING PORT WITH STERILE ALCOHOL PRIOR TO ACCESSING WITH A STERILE SYRINGE.**

PREPARATION PROTOCOL:

<p>STEP 1:</p>  <p>Attach the sterile filter needle onto the VACLOK 60mL syringe. Then draw 15mL of Heparin Anticoagulant (1000 units/mL)</p>	<p>STEP 2:</p>  <p>Then prime the bone marrow aspirating cannula by injecting 5mL of heparin through it</p>	<p>STEP 3:</p>  <p>Attach to the OUT port of the bone marrow filter. Fill to prime and then aspirate back into the syringe</p>	<p>STEP 4:</p>  <p>Then discard the residual heparin leaving 5mL in each VACLOK syringe</p>
--	---	---	---

CONCENTRATING PROTOCOL

<p>Step 1:</p>  <p>Slowly draw 55mL of bone marrow aspirate, filling each syringe to 60mL. Mix the BMA and heparin upon collection to prevent coagulation.</p>	<p>Step 2:</p>  <p>Start End</p> <p>For each BMA syringe filter the BMA in the direction indicated on the filter.</p>	<p>Step 3:</p>  <p>Inject 60mL of anticoagulated & filtered BMA into each Concentrating Device</p>	<p>Step 4:</p>  <p>Counterbalance with the same volume in each Concentrating Device.</p>  <p>Then place them directly opposite to each other in the centrifuge rotor buckets.</p>
--	---	--	--


STEP 5:

Sapphire Series
Centrifuge:
AbsolutePRP/BMC

Platinum Series
Centrifuge:
AbsolutePRP/BMC


Executive Series
Centrifuge:
10 minutes
4400 RPM

STEP 6:



Attach the 60mL syringe and aspirate the plasma until the bone marrow RBC interface reaches the 7mL volume marker, then stop aspirating

REPEAT FOR EACH CONCENTRATING DEVICE:



FOR EACH DEVICE

Attach the 12mL syringe and aspirate 7mL BMC
(Rotate the syringe to re-suspend the BMC buffycoat into the plasma)

COLLECT A TOTAL OF 14mL BMC