

# AbsolutePRP™ *Gold*



Hematocrit	Less than 1%
Avg Platelet Concentration	2.3
Avg Platelet Yield	92%
Avg Deliverable Platelets	337,808,333
Blood Draw Volume	20mL
Avg PRP Volume	8.2mL



## The AbsolutePRP™ Gold Concentrating System

The AbsolutePRP™ Gold Concentrating System is the resolute solution to high platelet yielding PRP from a low volume blood draw of 20mL. It's made of durable medical grade polycarbonate polymers that permits the clear and natural view of the PURE platelet cloud, that's devoid of red blood cells after centrifugation. Its unique design characteristics make it the fastest PRP system available, preparing PRP in one minute processing time. It is also the most affordable, high quality PRP that offers a therapeutic dose of growth factors, mediators and cytokines at the point of care. AbsolutePRP™ Gold is the best clinical option for low blood draw PRP applications.